REPROCESSING INSTRUCTIONS

"I FEEL
GOOD"

EMS 

EN: REPROCESSING INSTRUCTIONS

This IFU is only valid for North America.
This IFU is not valid for Australia. Please refer to FB-358/AU.
This IFU is not valid for China. Please refer to FB-358/CN.
Refer to FB-358 for the rest of the world.

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SUMMARY

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BEFORE USE

Please read these Operating Instructions carefully as they explain all the most important details and procedures. Please pay special attention to the safety precautions. Always keep these Instructions close at hand.

⚠ These Instructions are only applicable to the equipment they were delivered with.

To prevent injury to persons or damage to property, please observe the corresponding directives and symbols.

Translations in this manual are based on the English version. Thus the English text will be referred to in case of discrepancy.

DEVICES AND HANDPIECE CORDS


⚠ The following instructions have been validated for all reusable products¹ of EMS dental devices. It remains the responsibility of the processor to ensure that the processing, as actually performed using the equipment, materials and personnel in the processing facility, achieves the desired result. This requires validation and/or validation and routine monitoring of the process. Likewise any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse effects.

AIRFLOW® and PIEZON® devices

Handy

PIEZON® Kit on your unit

CLEAN

Wipe the surface with Clean the external surface of the device with a compatible EPA-registered intermediate-level surface disinfectant.

DISINFECT

Wipe the surface with 70% isopropyl alcohol or a disinfection wipe with proven efficacy.

⚠ Follow carefully the instructions provided by the disinfection solution manufacturer.

⚠ Do not use CaviWipes™. It damages EMS products.

INTRODUCTION

EMS recommends compliance with cleaning, packaging for sterilization and sterilization procedures according to ISO 17664.

⚠ Always report adverse events related to device reprocessing directly to EMS.

⚠ Products must be cleaned and, if applicable, sterilized prior to first use. Do not reprocess the products over the authorized number of sterilization cycles. To replace, refer to the "Service life" section of these Instructions.

⚠ Concentrations and contact times specified by the manufacturer of the cleaning agent must be observed.

⚠ Remember that sterilization cannot be achieved unless assembly components are cleaned first.

If parts of these instructions are not clear or seem to be inadequate, do not hesitate to contact/inform EMS.

⚠ Users shall also observe any legal requirements applicable in their country as well as the hygiene regulations in effect at the hospital or clinic.

This applies especially with regard to additional requirements for the inactivation of prions.

SERVICE LIFE

The products have been designed for a large number of sterilization cycles. The materials used in their manufacture were selected accordingly. However, with every renewed preparation for use, thermal and chemical stresses will result in ageing of the products.

⚠ Always replace products that present signs of wear or early degradation, regardless of the number of sterilization cycles left unused.

⚠ DO NOT expose the products to temperature exceeding 280°F (138°C).

⚠ EMS handpieces must not be lubricated or otherwise exposed to oil whether directly or indirectly.

Lifetime : see product's Instructions for Use.

SYMBOLS

⚠ Caution / Danger
Risk of property damage or injury

ⓘ Please note
Useful additional information and hints

⊘ Prohibited

🧤🧴 Wear personal protective equipment

📖 Refer to instruction manual

135°C 🔄 Sterilizable at up to 135°C in the autoclave

🧼 Suitable for thermal disinfection

⌚ Required time

Abbreviations used: min (minute) ; s (second)

PROCESS

PREPARE

⌚ 20 s - Immediately after use, rinse the lumen(s) line of the handpiece/instrument with water. Coarse soiling must be removed immediately after application.

Safely transport to the reprocessing area to avoid any damage to the products, contamination of the environment or of the persons involved in the reprocessing process.

⚠ ⌚ Cleaning is to be performed within 1 hour from use.

🧤🧴 Wear personal protective equipment according to the type of preparation.

! AIRFLOW® and PERIOFLOW® Handpieces:

Always carry out handpiece powder unclogging and check for both lumens (water and powder) clearance before proceeding. Use Easy Clean.

! PIEZON® Endochuck:

Remove any installed file.
Remove the screw and separate the small O-ring.

! PIEZON® Handpieces:

Remove any installed instrument.
Remove the nose cap and separate the light guide and O-ring.

! PIEZON® Instruments with CombiTorque®:

Separate aside the Instrument and the CombiTorque®.

Any product can be cleaned manually or automatically by washer or disinfector.

⚠ EMS handpieces must not be lubricated or otherwise exposed to oil whether directly or indirectly.

AUTOMATIC

🧴🌀 CLEAN AND DRY

Correctly place the product into a suitable rack, connect all lumens to the rinsing connectors and start the automated cleaning.

ⓘ Instructions of the manufacturer of the washer disinfector must also be observed.

See process in the following page.

MANUAL

🧴 CLEAN

The 2 following processes without ultrasound bath can be used with any EMS product:

Cleaning (without ultrasound bath) with Enzymax:

- Wipe the products using a lint-free cloth dampened with tap water to remove gross soil.
- Special attention need to be paid to cracks, crevices, seams and hard to reach areas.
- Prepare a bath of Enzymax cleaning solution following the manufacturer's recommendation of 1 oz/gallon of lukewarm tap water.
- Use a soft bristled brush to brush the article in the prepared cleaning solution until all visible soils have been removed.

- ⌚ 15 min Immerse the article in the prepared detergent solution. While immersed, ensure that all lumens are filled with the cleaning solution. Use a syringe, if needed. All surfaces must be moistened.
- ⌚ 15 s Flush each lumen with RO/DI water using a spray gun (water jet gun, with static water pressure of 2 bar).
- ⌚ 10 s Rinse the whole part under running RO/DI water.

- Fully dry the lumen and whole part until no more residues of water were present (visible or detectable).

AUTOMATIC

The following automated process¹ can be used:

- ⌚ 2 min Pre-washing with cold tap water².
Drain
- ⌚ 5 min Clean at 55°C with tap water² and 0.5%
cleaning solution³.
Drain
- ⌚ 3 min Rinse and neutralize with cold
demineralized water⁴.
Drain
- ⌚ 2 min Rinse with cold demineralized water⁴.
Drain

MANUAL

Cleaning (without ultrasound bath) with Enzol 2%:

- Brush carefully with a soft surface brush or cloth to remove any externally attached soiling.
- Remove the handpiece from the connector. Remove the gaskets, caps and other removable parts.

- ⌚ 10 min Soak all parts in the cleaning solution for, as a minimum, the time and the concentration specified by the manufacturer of the cleaning agent (ENZOL, 2 %, 10 min).

- Brush all lumens and surfaces with a suitable lumen / surface brush until all visible soiling is removed.
- Rinse through all product lumina (e.g. irrigation and aspiration connection) with ENZOL, 2 %, at least 5 times in the flow direction (no back rinsing), using a disposable syringe (min. volume 50 ml) applied to the nozzles of the product.
- Submerge all components completely in deionized water.
- Rinse through all product lumina (e.g. irrigation and aspiration connection) with deionized water at least 3 times in the flow direction (no back rinsing), using a disposable syringe (min. volume 50 ml) applied to the nozzles of the product.
- Repeat the cleaning process if the last rinse does not run clear, or if stains are still visible on the product.
- Dry all components at room temperature.

DO NOT use any ultrasound bath cleaning procedure with the handpieces: it may destroy the products.

The following process is intended **ONLY** for EMS PIEZON® Instruments and file holders:

Cleaning (with ultrasound bath) with Enzol 2%:

- Remove any externally attached soiling by brushing carefully with a soft surface brush or cloth.
- Remove the instruments from handpieces.

- ⌚ 10 min Soak all parts in the cleaning solution for, at a minimum, the time and the concentration specified by the manufacturer of the cleaning agent (ENZOL, 2 %, 10 min).

- Remove any externally attached soiling by brushing carefully with a soft surface brush or cloth.
- Place all parts in an ultrasound bath filled with ENZOL, 2 %.

- ⌚ 3 min (at least) Treat the parts with ultrasound

- Brush the parts in ENZOL, 2 % if any visible soiling is detected.
- Submerge all components completely in deionized water.
- Rinse through all product lumens with deionized water.
- Repeat the cleaning process if the last rinse does not run clear, or if stains are still visible on the product.
- Dry all components at room temperature.

DRY

Dry all products at room temperature using a clean lint free cloth.

INSPECT

! If stains are still visible on the product after cleaning, the entire cleaning procedure must be repeated. Products with visible damage, chip/flake loss, corrosion or deformation must be disposed of (no further use is allowed). Check also the integrity of O-rings and gaskets and replace if damaged or deformed.

! Verify the product to be fully dry. In case of detection of residues of water, remove these using an air pistol (clean compressed air). Fully dry the lumen and the whole product until no more residues of water are present (visible or detectable).

PACK

! Only previously cleaned products can be sterilized.

! Effective sterilization can take place only on fully dry products. Ensure each part (internal lumens and any surface) to be perfectly dry before reassembling and packing. Prior to sterilization, the products need to be reassembled for ready use and placed in suitable sterilization packaging. to their readiness of use and placed in a suitable sterilization packaging.


AIRFLOW® and PERIOFLOW® Handpieces:

No need for reassembly

PIEZON® Endochuck:

Reinstall the small O-ring (gasket) first, and then gently screw the retention bolt without tightening.

PIEZON® Handpieces:

Reinstall the O-ring first, then put the light guide into the nose cap and screw it on the handpiece.

PIEZON® Instruments with CombiTorque®:

Reinstall the Instrument on the CombiTorque®.

Any EMS part can be correctly packaged using an FDA-cleared sterilization pouch that has been validated for the specified sterilization cycle.

STERILIZE

! Sterilization must be performed immediately after cleaning.

DO NOT exceed

- the maximum number of sterilization cycles allowed.
- a sterilization temperature of 280°F (138°C) and a holding time of 20 min.

DO NOT use a class N sterilizer to reprocess PIEZON® handpieces.

DO NOT use hot-air sterilization and radio-sterilization procedures: they destroy the products.

Moist heat sterilization of products shall be performed according to ISO 17665 and under consideration of the respective country requirements.

! EMS recommends the use of biological indicators and chemical indicators that have been validated for use in a 132°C pre-vacuum steam sterilization cycle.

The prevacuum moist heat (steam) process can be used with any EMS product packaged in appropriate single or double pouches:

Parameters for the prevacuum moist heat cycle:

- 3 prevacuum phases
- Pressure of 3 bar¹
- Humidity of 100%
- Temperature of 270°F (132°C)
- ⌚ 4 min (minimum) Hold time (full cycle)
- ⌚ 20 min (minimum) Dry

! It is the duty of the user to ensure that the reprocessing processes, including resources, materials and personnel, are capable to reach the required results and maintained over the time: keeping actual the validation of the reprocessing processes is under the responsibility of the user.

STORE

! Store the sterilized products at a temperature of 41°F to 104°F (5°C to 40°C) in a:

- dry,
- clean
- and dust-free environment.

¹ For example: Miele Professional G 7836 CD with Miele Rack E429

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FROM THE VALLEY


TO THE TOP

► For more than 35 years EMS has developed high-end technologies and protocols in conjunction with the world's most advanced clinicians, enabling the patients to enjoy natural teeth and implants for much longer.

„Primum non nocere“ and minimal invasiveness are not only a must, they are our way of life.

MONT BLANC, 4807M

EMS RESEARCH AND MANUFACTURE CENTER
NYON, SWITZERLAND

ONLY USE ORIGINAL EMS PRODUCTS

► EMS powders are much less abrasive than traditional mechanical debridement techniques such as rubber cups and handscaling.
We recommend to use AIRFLOW® and PERIOFLOW® first, followed by PIEZON®. GBT: the minimally invasive way for Professional Mechanical Plaque Removal (PMPR).

► As part of the GBT disclosing phase, coloration also minimizes AIRFLOW® abrasion. No more color signifies that all biofilm has been removed: the objective of AIRFLOW® is achieved, to preserve the integrity of dental tissues.

⚠ EMS devices have been designed and tested for exclusive use with EMS prophylaxis powders.

! ONLY USE EMS POWDERS

So-called "EMS compatible" powders in the marketplace may cause damage to oral tissues and the device. See below a damaged Airflow nozzle after use with a "compatible" powder. If EMS products are damaged by using "compatible" powders the product Warranty will be lost.

Damaged handpiece

ONLY USE THE SWISS ORIGINAL EMS INSTRUMENTS!

The EMS Instrument, the Handpiece and the NO PAIN electronic module have to vibrate in harmony like in a trilogy.
Using so-called "compatible" and copy tips from third parties can damage your patients' teeth as well as EMS handpieces.
They were developed and made for each other by EMS.
One PS Instrument is designed for up to 500 treatments, thus reducing the cost of one treatment to only a few cents.
Therefore, there is absolutely no reason to "save" money by using so-called "compatible" tips. Moreover, if the Piezon Handpiece gets damaged (the thread) it will lose its warranty so at the end of the day, "compatible" tips will always cost you more.

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AU

SYDNEY, AUSTRALIA
E.M.S. ELECTRO MEDICAL SYSTEMS OCEANIA
Unit 1, Building 4 – Sydney Corporate Park
160 Bourke Road
AU - Alexandria, NSW 2015
T. +61 405 095 867
E-mail: info@ems-oceania.com

JP

TOKYO, JAPAN
E.M.S. JAPAN KK
Eminent Ningyoucho Bld 6F.,
9-18, Kobuna-cho, Nihombashi
JP - Chuo-ku, Tokyo, 103-0024
T. +81 (0) 3-6810-9351
F. +81 (0) 3-6810-9352
E-mail: emsjapan@ems-ch.com

CN

SHANGHAI, CHINA
医迈斯电子医疗系统贸易 (上海) 有限公司
E.M.S. ELECTRO MEDICAL SYSTEMS TRADING (SHANGHAI) CO., LTD.
24A, Jin Sui Mansion,
N°379 Pudong Nan Rd.
CN - Shanghai 200120
T. +86 21 3363 2323
E-mail: emschina@ems-ch.com

KR

SEOUL, SOUTH KOREA
EMS ELECTRO MEDICAL SYSTEMS KOREA LTD..
13, Beobwon-ro 8-gil, Songpa-gu
T. +82 70 7807 2875
E-mail: emskorea@ems-ch.com

EMS 

WWW.EMS-DENTAL.COM